

KERN 315-312-407

Safety Set 100 g (M1); 2000 g (F2)

KERN



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	Safety Set
Product family	315

Measuring System

Test weight - nominal value	100 g & 2 kg
Weight OIML class	F2 M1

Construction

Number of weights	2
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	255×135×140 mm
Net weight	0,01 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	1 kg

Services (optional)

Article number for DAkkS calibration	962-400-600
--------------------------------------	-------------