

KERN 316-01

Individual weight, OIML Class E2, 1 g, Compact form, stainless steel polished (OIML)



Services (optional)

| | |
|--|---------|
| Article number for DAkkS calibration | 962-331 |
| Article number for guaranteed plus-deviation | 921-331 |

Pictograms

STANDARD



Category

| | |
|------------------|-------------------|
| Brand | KERN |
| Product category | Individual weight |
| Product group | OIML E2 |
| Product family | 316 |

Measuring System

| | |
|---|---------|
| Test weight - nominal value | 1 g |
| Test weight - volume (cm ³) | 0,125 |
| Test weight - volume uncertainty (cm ³) | 0,0016 |
| Permissible tolerance according to OIML (+/-) | 0,03 mg |
| Test weight - Density [Min] (OIML) | 5300 |
| Test weight - Density [Max] (OIML) | 16000 |
| Weight OIML class | E2 |

Construction

| | |
|------------------|---------------------------------|
| Shape weight | Compact form |
| Dimensions (Ø×H) | 6×6 mm |
| Material | stainless steel polished (OIML) |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (W×D×H) | 6×6×6 mm |
| Net weight | 0,001 kg |
| Shipping method | Parcel service |
| Net weight approx. | 0,05 kg |
| Gross weight approx. | 0,05 kg |
| Shipping weight | 0,1 kg |