

KERN 316-09

KERN

Individual weight, OIML Class E2, 500 g, Compact form, stainless steel polished (OIML)



Services (optional)

Article number for DAkkS calibration	962-339
Article number for guaranteed plus-deviation	921-339

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	316

Measuring System

Test weight - nominal value	500 g
Test weight - volume (cm ³)	62,5
Test weight - volume uncertainty (cm ³)	0,7813
Permissible tolerance according to OIML (+/-)	0,8 mg
Test weight - Density [Min] (OIML)	7810
Test weight - Density [Max] (OIML)	8210
Weight OIML class	E2

Construction

Shape weight	Compact form
Dimensions (Ø×H)	38×58,5 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	38×38×58 mm
Net weight	0,5 kg
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,50 kg
Shipping weight	0,5 kg