

# KERN 316-12

**KERN**

Individual weight, OIML Class E2, 2 kg, Compact form, stainless steel polished (OIML)



## Services (optional)

Article number for DAkkS calibration	962-342
Article number for guaranteed plus-deviation	921-342

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	316

## Measuring System

Test weight - nominal value	2 kg
Test weight - volume (cm <sup>3</sup> )	250
Test weight - volume uncertainty (cm <sup>3</sup> )	3,125
Permissible tolerance according to OIML (+/-)	3 mg
Test weight - Density [Min] (OIML)	7810
Test weight - Density [Max] (OIML)	8210
Weight OIML class	E2

## Construction

Shape weight	Compact form
Dimensions (Ø×H)	60×94 mm
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	60×60×94 mm
Net weight	2 kg
Shipping method	Parcel service
Net weight approx.	2,0 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg