

KERN 317-020-600

KERN

Aluminium protective box for individual weights (Cylindrical/compact, 2 g)



Category

| | |
|------------------|--------------------------|
| Brand | KERN |
| Product category | Cases & Boxes |
| Product group | Aluminium protective box |

Measuring System

| | |
|-----------------------------------------|---------------------|
| Case suitable for OIML classes | E1 - M3 |
| Case suitable for nominal weight values | 2 g |
| Suitable weight shape for case | Cylindrical/compact |

Construction

| | |
|-------------------|-------------|
| Dimension (W×D×H) | 70×77×73 mm |
| Material | aluminium |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (W×D×H) | 70×77×73 mm |
| Net weight | 0,18 kg |
| Shipping method | Parcel service |
| Net weight approx. | 0,20 kg |
| Gross weight approx. | 0,20 kg |
| Shipping weight | 0,18 kg |

Pictograms

STANDARD

