

KERN 317-060-400

KERN

Plastic box for individual weights (Cylindrical/compact, 50 g)



Category

| | |
|------------------|---------------|
| Brand | KERN |
| Product category | Cases & Boxes |
| Product group | Plastic box |

Measuring System

| | |
|---|---------------------|
| Case suitable for OIML classes | E1, E2 |
| Case suitable for nominal weight values | 50 g |
| Suitable weight shape for case | Cylindrical/compact |

Construction

| | |
|-----------------|----------|
| Material | plastic |
| Dimension (Ø×H) | 31×56 mm |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (W×D×H) | 31×31×56 mm |
| Net weight | 0,02 kg |
| Shipping method | Parcel service |
| Net weight approx. | 0,05 kg |
| Gross weight approx. | 0,05 kg |
| Shipping weight | 0,1 kg |

Pictograms

STANDARD

