

KERN 317-11

Individual weight, OIML Class E2, 1 kg, Knob, stainless steel polished (OIML)



Services (optional)

Article number for conformity assessment	952-341
Article number for DAkkS calibration	962-341
Article number for guaranteed plus-deviation	921-341

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	317

Measuring System

Test weight - nominal value	1 kg
Test weight - volume (cm ³)	125
Test weight - volume uncertainty (cm ³)	1,5625
Permissible tolerance according to OIML (+/-)	1,6 mg
Test weight - Density [Min] (OIML)	7810
Test weight - Density [Max] (OIML)	8210
Weight OIML class	E2

Construction

Shape weight	Knob
Dimensions (Ø×H)	48×82 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	48×48×82 mm
Net weight	1 kg
Shipping method	Parcel service
Net weight approx.	1,00 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg