

KERN 317-12

KERN

Individual weight, OIML Class E2, 2 kg, Knob, stainless steel polished (OIML)



Services (optional)

Article number for conformity assessment	952-342
Article number for DAkkS calibration	962-342
Article number for guaranteed plus-deviation	921-342

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	317

Measuring System

Test weight - nominal value	2 kg
Test weight - volume (cm ³)	250
Test weight - volume uncertainty (cm ³)	3,125
Permissible tolerance according to OIML (+/-)	3 mg
Test weight - Density [Min] (OIML)	7810
Test weight - Density [Max] (OIML)	8210
Weight OIML class	E2

Construction

Shape weight	Knob
Dimensions (Ø×H)	60×100 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	60×60×100 mm
Net weight	2 kg
Shipping method	Parcel service
Net weight approx.	2,0 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg