

KERN 317-13

Individual weight, OIML Class E2, 5 kg, Knob, stainless steel polished (OIML)



Services (optional)

Article number for conformity assessment	952-343
Article number for DAkkS calibration	962-343
Article number for guaranteed plus-deviation	921-343

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	317

Measuring System

Test weight - nominal value	5 kg
Test weight - volume (cm ³)	625
Test weight - volume uncertainty (cm ³)	7,8125
Permissible tolerance according to OIML (+/-)	8 mg
Test weight - Density [Min] (OIML)	7810
Test weight - Density [Max] (OIML)	8210
Weight OIML class	E2

Construction

Shape weight	Knob
Dimensions (Ø×H)	80×140 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	80×80×140 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg