

KERN 317-130-100



Wooden box for individual weights (Cylindrical/compact, 5 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	E1+E2+F1
Case suitable for nominal weight values	5 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	106×106×178 mm
Material	beech

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	106×106×178 mm
Net weight	0,57 kg
Shipping method	Parcel service
Net weight approx.	0,60 kg
Gross weight approx.	0,60 kg
Shipping weight	0,58825 kg

Pictograms

STANDARD

