

# KERN 317-15

Individual weight, OIML Class E2, 20 kg, Knob, stainless steel polished (OIML)



## Services (optional)

|  |         |
|--|---------|
| Article number for conformity assessment     | 952-345 |
| Article number for DAkkS calibration         | 962-345 |
| Article number for guaranteed plus-deviation | 921-345 |

## Pictograms

### STANDARD



## Category

|                  |                   |
|------------------|-------------------|
| Brand            | KERN              |
| Product category | Individual weight |
| Product group    | OIML E2           |
| Product family   | 317               |

## Measuring System

|   |       |
|---|-------|
| Test weight - nominal value                         | 20 kg |
| Test weight - volume (cm <sup>3</sup> )             | 2500  |
| Test weight - volume uncertainty (cm <sup>3</sup> ) | 31,25 |
| Permissible tolerance according to OIML (+/-)       | 30 mg |
| Test weight - Density [Min] (OIML)                  | 7810  |
| Test weight - Density [Max] (OIML)                  | 8210  |
| Weight OIML class                                   | E2    |

## Construction

|                  |                                 |
|------------------|---------------------------------|
| Shape weight     | Knob                            |
| Dimensions (Ø×H) | 129×227 mm                      |
| Material         | stainless steel polished (OIML) |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Readability force [d] (N)    | 1 d            |
| Dimensions packaging (W×D×H) | 129×129×227 mm |
| Net weight                   | 20 kg          |
| Shipping method              | Parcel service |
| Net weight approx.           | 20 kg          |
| Gross weight approx.         | 20 kg          |
| Shipping weight              | 20 kg          |