

KERN 318-226

KERN

Set of milligram weights, E2, 1 mg - 500 mg, Platelet, stainless steel polished (OIML), in aluminium case



Article number for guaranteed plus-deviation 921-350

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	318

Measuring System

Test weight - nominal value	1 mg - 500 mg
Weight OIML class	E2

Construction

Shape weight	Platelet
Dimension (W×D×H)	165×158×37 mm
Number of weights	13
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	165×158×37 mm
Net weight	0,31 kg
Shipping method	Parcel service
Net weight approx.	0,35 kg
Gross weight approx.	0,40 kg
Shipping weight	0,39 kg

Services (optional)

Article number for conformity assessment	952-350
Article number for DAkkS calibration	962-350