

# KERN 318-226

**KERN**

Set of milligram weights, E2, 1 mg - 500 mg, Platelet, stainless steel polished (OIML), in aluminium case



Article number for DAkkS calibration 962-350

Article number for guaranteed plus-deviation 921-350

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	318

## Measuring System

Test weight - nominal value	1 mg - 500 mg
OIML-compliant	✓
Weight OIML class	E2

## Construction

Shape weight	Platelet
Dimension (W×D×H)	165×158×37 mm
Number of weights	13
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	165×158×37 mm
Net weight	0,31 kg
Shipping method	Parcel service
Net weight approx.	0,35 kg
Gross weight approx.	0,40 kg
Shipping weight	0,39 kg

## Services (optional)

Article number for conformity assessment	952-350
--	---------