

# KERN 318-37

Milligram weights, E2, 100 mg, Wire, stainless steel polished (OIML)



## Services (optional)

Article number for conformity assessment	952-357
Article number for DAkkS calibration	962-357
Article number for guaranteed plus-deviation	921-357

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	318

## Measuring System

Test weight - nominal value	100 mg
Test weight - volume (cm <sup>3</sup> )	0,0125
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0003
Permissible tolerance according to OIML (+/-)	0,016 mg
Weight OIML class	E2

## Construction

Shape weight	Wire
Dimensions weight (W×D×H)	9×9×3 mm
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	9×9×3 mm
Net weight	0,0001 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg