

# KERN 323-024

# KERN

Weight set, F1, 1 mg - 50 g, Knob, stainless steel polished (OIML), in plastic case



Article number for guaranteed plus-deviation 921-401

Pictograms

STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

## Measuring System

Test weight - nominal value	1 mg - 50 g
Weight OIML class	F1

## Construction

Shape weight	Knob
Dimension (W×D×H)	250×210×75 mm
Number of weights	20
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Net weight	0,5 kg
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,50 kg
Shipping weight	0,8 kg

## Services (optional)

Article number for conformity assessment	952-401
Article number for DAkkS calibration	962-401