

KERN 323-03

KERN

Weight set, F1, 1 mg - 100 g, Knob, stainless steel polished (OIML), in wooden case



Article number for guaranteed plus-deviation 921-402

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

Measuring System

Test weight - nominal value	1 mg - 100 g
Weight OIML class	F1

Construction

Shape weight	Knob
Dimension (WxDxH)	152x152x66 mm
Number of weights	21
Case of weight set	in wooden case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	152x152x66 mm
Net weight	0,88 kg
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	0,90 kg
Shipping weight	0,88 kg

Services (optional)

Article number for conformity assessment	952-402
Article number for DAkkS calibration	962-402