

KERN 323-036

KERN

Weight set, F1, 1 mg - 100 g, Knob, stainless steel polished (OIML), in aluminium case



Article number for guaranteed plus-deviation 921-402

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

Measuring System

Test weight - nominal value	1 mg - 100 g
Weight OIML class	F1

Construction

Shape weight	Knob
Dimension (WxDxH)	165x125x75 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	165x125x75 mm
Net weight	0,73 kg
Shipping method	Parcel service
Net weight approx.	0,75 kg
Gross weight approx.	0,90 kg
Shipping weight	0,9 kg

Services (optional)

Article number for conformity assessment	952-402
Article number for DAkkS calibration	962-402