

# KERN 323-04

# KERN

Weight set, F1, 1 mg - 200 g, Knob, stainless steel polished (OIML), in wooden case



Article number for guaranteed plus-deviation 921-403

Pictograms

STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

## Measuring System

Test weight - nominal value	1 mg - 200 g
Weight OIML class	F1

## Construction

Shape weight	Knob
Dimension (W×D×H)	152×152×66 mm
Number of weights	23
Case of weight set	in wooden case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	152×152×66 mm
Net weight	1,29 kg
Shipping method	Parcel service
Net weight approx.	1,4 kg
Gross weight approx.	1,4 kg
Shipping weight	1,29 kg

## Services (optional)

Article number for conformity assessment	952-403
Article number for DAkkS calibration	962-403