

KERN 323-06

KERN

Weight set, F1, 1 mg - 1 kg, Knob, stainless steel polished (OIML), in wooden case



Article number for guaranteed plus-deviation 921-405

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

Measuring System

Test weight - nominal value	1 mg - 1 kg
Weight OIML class	F1

Construction

Shape weight	Knob
Dimension (WxDxH)	230x180x110 mm
Number of weights	25
Case of weight set	in wooden case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	230x180x110 mm
Net weight	4,21 kg
Shipping method	Parcel service
Net weight approx.	4,4 kg
Gross weight approx.	4,4 kg
Shipping weight	4,21 kg

Services (optional)

Article number for conformity assessment	952-405
Article number for DAkkS calibration	962-405