

KERN 323-064

KERN

Weight set, F1, 1 mg - 1 kg, Knob, stainless steel polished (OIML), in plastic case



Article number for guaranteed plus-deviation 921-405

Pictograms

STANDARD



Category

| | |
|------------------|------------|
| Brand | KERN |
| Product category | Weight set |
| Product group | OIML F1 |
| Product family | 323 |

Measuring System

| | |
|-----------------------------|-------------|
| Test weight - nominal value | 1 mg - 1 kg |
| Weight OIML class | F1 |

Construction

| | |
|--------------------|---------------------------------|
| Shape weight | Knob |
| Dimension (W×D×H) | 350×320×180 mm |
| Number of weights | 25 |
| Case of weight set | in plastic case |
| Material | stainless steel polished (OIML) |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (W×D×H) | 350×320×180 mm |
| Net weight | 3,7 kg |
| Shipping method | Parcel service |
| Net weight approx. | 3,8 kg |
| Gross weight approx. | 3,8 kg |
| Shipping weight | 4 kg |

Services (optional)

| | |
|--|---------|
| Article number for conformity assessment | 952-405 |
| Article number for DAkkS calibration | 962-405 |