

# KERN 323-066

# KERN

Weight set, F1, 1 mg - 1 kg, Knob, stainless steel polished (OIML), in aluminium case



Article number for guaranteed plus-deviation 921-405

Pictograms

STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

## Measuring System

Test weight - nominal value	1 mg - 1 kg
Weight OIML class	F1

## Construction

Shape weight	Knob
Dimension (W×D×H)	230×140×110 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	230×140×110 mm
Net weight	2,78 kg
Shipping method	Parcel service
Net weight approx.	2,8 kg
Gross weight approx.	3,6 kg
Shipping weight	3,5 kg

## Services (optional)

Article number for conformity assessment	952-405
Article number for DAKkS calibration	962-405