

# KERN 323-074

**KERN**

Weight set, F1, 1 mg - 2 kg, Knob, stainless steel polished (OIML), in plastic case



Article number for guaranteed plus-deviation 921-406

Pictograms

STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

## Measuring System

Test weight - nominal value	1 mg - 2 kg
Weight OIML class	F1

## Construction

Shape weight	Knob
Dimension (W×D×H)	350×320×180 mm
Number of weights	27
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	350×320×180 mm
Net weight	7,69 kg
Shipping method	Parcel service
Net weight approx.	8 kg
Gross weight approx.	8 kg
Shipping weight	7,69 kg

## Services (optional)

Article number for conformity assessment	952-406
Article number for DAkkS calibration	962-406