

KERN 323-084

KERN

Weight set, F1, 1 mg - 5 kg, Knob, stainless steel polished (OIML), in plastic case



Article number for guaranteed plus-deviation 921-407

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

Measuring System

Test weight - nominal value	1 mg - 5 kg
Weight OIML class	F1

Construction

Shape weight	Knob
Dimension (W×D×H)	350×320×180 mm
Number of weights	28
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	350×320×180 mm
Net weight	12,69 kg
Shipping method	Parcel service
Net weight approx.	13 kg
Gross weight approx.	13 kg
Shipping weight	12,69 kg

Services (optional)

Article number for conformity assessment	952-407
Article number for DAkkS calibration	962-407