

KERN 324-024

KERN

Weight set, F1, 1 g - 50 g, Knob, stainless steel polished (OIML), in plastic case



Article number for guaranteed plus-deviation 921-415

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	324

Measuring System

Test weight - nominal value	1 g - 50 g
Weight OIML class	F1

Construction

Shape weight	Knob
Dimension (W×D×H)	250×210×75 mm
Number of weights	8
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Net weight	0,5 kg
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,50 kg
Shipping weight	0,8 kg

Services (optional)

Article number for conformity assessment	952-415
Article number for DAkkS calibration	962-415