

# KERN 324-04

**KERN**

Weight set, F1, 1 g - 200 g, Knob, stainless steel polished (OIML), in wooden case



Article number for DAkkS calibration 962-417

Article number for guaranteed plus-deviation 921-417

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	324

## Measuring System

Test weight - nominal value	1 g - 200 g
OIML-compliant	✓
Weight OIML class	F1

## Construction

Shape weight	Knob
Dimension (W×D×H)	152×152×66 mm
Number of weights	11
Case of weight set	in wooden case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	152×152×66 mm
Net weight	1,29 kg
Shipping method	Parcel service
Net weight approx.	1,4 kg
Gross weight approx.	1,4 kg
Shipping weight	1,29 kg

## Services (optional)

Article number for conformity assessment	952-417
--	---------