

KERN 324-076

KERN

Weight set, F1, 1 g - 2 kg, Knob, stainless steel polished (OIML), in aluminium case



Article number for guaranteed plus-deviation 921-420

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	324

Measuring System

Test weight - nominal value	1 g - 2 kg
Weight OIML class	F1

Construction

Shape weight	Knob
Dimension (WxDxH)	278x180x138 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	278x180x138 mm
Net weight	5,15 kg
Shipping method	Parcel service
Net weight approx.	6 kg
Gross weight approx.	6 kg
Shipping weight	5,55 kg

Services (optional)

Article number for conformity assessment	952-420
Article number for DAKkS calibration	962-420