

KERN 324-086

Weight set, F1, 1 g - 5 kg, Knob, stainless steel polished (OIML), in aluminium case



Article number for guaranteed plus-deviation 921-421

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	324

Measuring System

Test weight - nominal value	1 g - 5 kg
Weight OIML class	F1

Construction

Shape weight	Knob
Dimension (W×D×H)	460×260×200 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	460×260×200 mm
Net weight	12,6 kg
Shipping method	Parcel service
Net weight approx.	13 kg
Gross weight approx.	18 kg
Shipping weight	18 kg

Services (optional)

Article number for conformity assessment	952-421
Article number for DAKkS calibration	962-421