

KERN 325-026

KERN

Weight set, F1, 1 mg - 50 g, Eco Shape, stainless steel polished (OIML), in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	325

Measuring System

Test weight - nominal value	1 mg - 50 g
Weight OIML class	F1

Construction

Shape weight	Eco Shape
Dimension (W×D×H)	165×120×75 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	165×120×75 mm
Net weight	0,54 kg
Shipping method	Parcel service
Net weight approx.	0,55 kg
Gross weight approx.	0,80 kg
Shipping weight	0,8 kg

Services (optional)

Article number for DAkkS calibration	962-401
Article number for guaranteed plus-deviation	921-401