

KERN 325-034

KERN

Weight set, F1, 1 mg - 100 g, Eco Shape, stainless steel polished (OIML), in plastic case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	325

Measuring System

Test weight - nominal value	1 mg - 100 g
Weight OIML class	F1

Construction

Shape weight	Eco Shape
Dimension (W×D×H)	250×210×75 mm
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Net weight	0,53 kg
Shipping method	Parcel service
Net weight approx.	0,55 kg
Gross weight approx.	0,60 kg
Shipping weight	0,8 kg

Services (optional)

Article number for DAkkS calibration	962-402
Article number for guaranteed plus-deviation	921-402