

# KERN 325-046

**KERN**

Weight set, F1, 1 mg - 200 g, Eco Shape, stainless steel polished (OIML), in aluminium case



## Pictograms

### STANDARD



## Category

|                  |            |
|------------------|------------|
| Brand            | KERN       |
| Product category | Weight set |
| Product group    | OIML F1    |
| Product family   | 325        |

## Measuring System

|                             |              |
|-----------------------------|--------------|
| Test weight - nominal value | 1 mg - 200 g |
| Weight OIML class           | F1           |

## Construction

|                    |                                 |
|--------------------|---------------------------------|
| Shape weight       | Eco Shape                       |
| Dimension (W×D×H)  | 190×130×80 mm                   |
| Case of weight set | in aluminium case               |
| Material           | stainless steel polished (OIML) |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Readability force [d] (N)    | 1 d            |
| Dimensions packaging (W×D×H) | 190×130×80 mm  |
| Net weight                   | 0,59 kg        |
| Shipping method              | Parcel service |
| Net weight approx.           | 0,60 kg        |
| Gross weight approx.         | 1,2 kg         |
| Shipping weight              | 1,1 kg         |

## Services (optional)

|  |         |
|--|---------|
| Article number for DAkkS calibration         | 962-403 |
| Article number for guaranteed plus-deviation | 921-403 |