

KERN 325-066

KERN

Weight set, F1, 1 mg - 1 kg, Eco Shape, stainless steel polished (OIML), in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	325

Measuring System

Test weight - nominal value	1 mg - 1 kg
Weight OIML class	F1

Construction

Shape weight	Eco Shape
Dimension (W×D×H)	230×140×110 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	230×140×110 mm
Net weight	0,89 kg
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	3,6 kg
Shipping weight	3,5 kg

Services (optional)

Article number for DAkkS calibration	962-405
Article number for guaranteed plus-deviation	921-405