

# KERN 325-076

**KERN**

Weight set, F1, 1 mg - 2 kg, Eco Shape, stainless steel polished (OIML), in aluminium case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	325

## Measuring System

Test weight - nominal value	1 mg - 2 kg
Weight OIML class	F1

## Construction

Shape weight	Eco Shape
Dimension (W×D×H)	278×180×138 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	278×180×138 mm
Net weight	1,34 kg
Shipping method	Parcel service
Net weight approx.	1,4 kg
Gross weight approx.	6 kg
Shipping weight	5,55 kg

## Services (optional)

Article number for DAkkS calibration	962-406
Article number for guaranteed plus-deviation	921-406