

KERN 325-092

KERN

Weight set, F1, 1 mg - 10 kg, Eco Shape, stainless steel polished (OIML), in wooden case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	325

Measuring System

Test weight - nominal value	1 mg - 10 kg
Weight OIML class	F1

Construction

Shape weight	Eco Shape
Dimension (W×D×H)	540×200×240 mm
Case of weight set	in wooden case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	540×200×240 mm
Net weight	32 kg
Shipping method	Freight
Net weight approx.	32 kg
Gross weight approx.	34 kg
Shipping weight	33 kg

Services (optional)

Article number for DAkkS calibration	962-408
Article number for guaranteed plus-deviation	921-408