

# KERN 326-034

**KERN**

Weight set, F1, 1 g - 100 g, Eco Shape, stainless steel polished (OIML), in plastic case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	326

## Measuring System

Test weight - nominal value	1 g - 100 g
Weight OIML class	F1

## Construction

Shape weight	Eco Shape
Dimension (W×D×H)	250×210×75 mm
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Net weight	0,36 kg
Shipping method	Parcel service
Net weight approx.	0,40 kg
Gross weight approx.	0,60 kg
Shipping weight	0,8 kg

## Services (optional)

Article number for DAkkS calibration	962-416
Article number for guaranteed plus-deviation	921-416