

# KERN 326-03

# KERN

Individual weight, OIML Class F1, 5 g, Eco Shape, stainless steel polished (OIML)

Article number for guaranteed plus-deviation 921-433

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

## Measuring System

Test weight - nominal value	5 g
Permissible tolerance according to OIML (+/-)	0,16 mg
Test weight - Density [Min] (OIML)	5300
Test weight - Density [Max] (OIML)	16000
Weight OIML class	F1

## Construction

Shape weight	Eco Shape
Dimensions (Ø×H)	8×13,2 mm
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	8×8×13 mm
Net weight	0,005 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Services (optional)

Article number for DAkkS calibration	962-433
--------------------------------------	---------