

KERN 326-046

KERN

Weight set, F1, 1 g - 200 g, Eco Shape, stainless steel polished (OIML), in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	326

Measuring System

Test weight - nominal value	1 g - 200 g
Weight OIML class	F1

Construction

Shape weight	Eco Shape
Dimension (W×D×H)	190×130×80 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	190×130×80 mm
Net weight	0,51 kg
Shipping method	Parcel service
Net weight approx.	0,55 kg
Gross weight approx.	1,2 kg
Shipping weight	1,1 kg

Services (optional)

Article number for DAkkS calibration	962-417
Article number for guaranteed plus-deviation	921-417