

KERN 326-06

KERN

Individual weight, OIML Class F1, 50 g, Eco Shape, stainless steel polished (OIML)



Article number for guaranteed plus-deviation 921-436

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

Measuring System

Test weight - nominal value	50 g
Permissible tolerance according to OIML (+/-)	0,3 mg
Test weight - Density [Min] (OIML)	7270
Test weight - Density [Max] (OIML)	8890
Weight OIML class	F1

Construction

Shape weight	Eco Shape
Dimensions (Ø×H)	18×29 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	18×18×29 mm
Net weight	0,05 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

Services (optional)

Article number for DAkkS calibration	962-436
--------------------------------------	---------