

# KERN 326-07

# KERN

Individual weight, OIML Class F1, 100 g, Eco Shape, stainless steel polished (OIML)



Article number for guaranteed plus-deviation 921-437

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

## Measuring System

Test weight - nominal value	100 g
Permissible tolerance according to OIML (+/-)	0,5 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

## Construction

Shape weight	Eco Shape
Dimensions (Ø×H)	22×38 mm
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	22×22×38 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

## Services (optional)

Article number for DAkkS calibration	962-437
--------------------------------------	---------