

KERN 326-08

KERN

Individual weight, OIML Class F1, 200 g, Eco Shape, stainless steel polished (OIML)



Article number for guaranteed plus-deviation 921-438

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

Measuring System

Test weight - nominal value	200 g
Permissible tolerance according to OIML (+/-)	1 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

Construction

Shape weight	Eco Shape
Dimensions (Ø×H)	28×47 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	28×28×47 mm
Net weight	0,2 kg
Shipping method	Parcel service
Net weight approx.	0,20 kg
Gross weight approx.	0,20 kg
Shipping weight	0,2 kg

Services (optional)

Article number for DAkkS calibration	962-438
--------------------------------------	---------