

# KERN 326-092

# KERN

Weight set, F1, 1 g - 10 kg, Eco Shape, stainless steel polished (OIML), in wooden case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	326

## Measuring System

Test weight - nominal value	1 g - 10 kg
Weight OIML class	F1

## Construction

Shape weight	Eco Shape
Dimension (W×D×H)	540×200×240 mm
Case of weight set	in wooden case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	540×200×240 mm
Net weight	30,61 kg
Shipping method	Freight
Net weight approx.	32 kg
Gross weight approx.	34 kg
Shipping weight	32,5 kg

## Services (optional)

Article number for DAkkS calibration	962-422
Article number for guaranteed plus-deviation	921-422