

# KERN 326-09

# KERN

Individual weight, OIML Class F1, 500 g, Eco Shape, stainless steel polished (OIML)



Article number for guaranteed plus-deviation 921-439

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

## Measuring System

Test weight - nominal value	500 g
Permissible tolerance according to OIML (+/-)	2,5 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

## Construction

Shape weight	Eco Shape
Dimensions (Ø×H)	38×64 mm
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	38×38×64 mm
Net weight	0,5 kg
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,50 kg
Shipping weight	0,5 kg

## Services (optional)

Article number for DAkkS calibration	962-439
--------------------------------------	---------