

KERN 326-12

Individual weight, OIML Class F1, 2 kg, Eco Shape, stainless steel polished (OIML)



Article number for guaranteed plus-deviation 921-442

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

Measuring System

Test weight - nominal value	2 kg
Permissible tolerance according to OIML (+/-)	10 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

Construction

Shape weight	Eco Shape
Dimensions (Ø×H)	60×102 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	60×60×102 mm
Net weight	2 kg
Shipping method	Parcel service
Net weight approx.	2,0 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg

Services (optional)

Article number for DAkkS calibration	962-442
--------------------------------------	---------