

KERN 326-13

Individual weight, OIML Class F1, 5 kg, Eco Shape, stainless steel polished (OIML)



Article number for guaranteed plus-deviation 921-443

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

Measuring System

Test weight - nominal value	5 kg
Permissible tolerance according to OIML (+/-)	25 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

Construction

Shape weight	Eco Shape
Dimensions (Ø×H)	80×144 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	80×80×144 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg

Services (optional)

Article number for DAkkS calibration	962-443
--------------------------------------	---------