

KERN 326-14

KERN

Individual weight, OIML Class F1, 10 kg, Eco Shape, stainless steel polished (OIML)



Article number for guaranteed plus-deviation 921-444

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

Measuring System

Test weight - nominal value	10 kg
Permissible tolerance according to OIML (+/-)	50 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

Construction

Shape weight	Eco Shape
Dimensions (Ø×H)	100×183 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	100×100×183 mm
Net weight	10 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg

Services (optional)

Article number for DAkkS calibration	962-444
--------------------------------------	---------