

KERN 326-39

KERN

Individual weight, OIML Class F1, 50 kg, Block, stainless steel polished (OIML)



Services (optional)

Article number for DAkkS calibration 962-446

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

Measuring System

Test weight - nominal value	50 kg
Test weight - volume uncertainty (cm ³)	78,125
Permissible tolerance according to OIML (+/-)	250 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

Construction

Shape weight	Block
Dimensions weight (W×D×H)	295×158×191 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	295×158×191 mm
Net weight	50 kg
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg