

KERN 327-03

KERN

Individual weight, OIML Class F1, 5 g, Knob, stainless steel polished (OIML)



Services (optional)

Article number for conformity assessment	952-433
Article number for DAkkS calibration	962-433
Article number for guaranteed plus-deviation	921-433

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

Measuring System

Test weight - nominal value	5 g
Test weight - volume (cm ³)	0,625
Test weight - volume uncertainty (cm ³)	0,0078
Permissible tolerance according to OIML (+/-)	0,16 mg
Test weight - Density [Min] (OIML)	5300
Test weight - Density [Max] (OIML)	16000
Weight OIML class	F1

Construction

Shape weight	Knob
Dimensions (Ø×H)	8×16,5 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	8×8×16 mm
Net weight	0,005 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg