

KERN 327-04

KERN

Individual weight, OIML Class F1, 10 g, Knob, stainless steel polished (OIML)



Services (optional)

Article number for conformity assessment	952-434
Article number for DAkkS calibration	962-434
Article number for guaranteed plus-deviation	921-434

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

Measuring System

Test weight - nominal value	10 g
Test weight - volume (cm ³)	1,25
Test weight - volume uncertainty (cm ³)	0,0156
Permissible tolerance according to OIML (+/-)	0,2 mg
Test weight - Density [Min] (OIML)	6000
Test weight - Density [Max] (OIML)	12000
Weight OIML class	F1

Construction

Shape weight	Knob
Dimensions (Ø×H)	10×19 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	10×10×19 mm
Net weight	0,01 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg