

KERN 327-07

KERN

Individual weight, OIML Class F1, 100 g, Knob, stainless steel polished (OIML)



Services (optional)

Article number for conformity assessment	952-437
Article number for DAkkS calibration	962-437
Article number for guaranteed plus-deviation	921-437

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

Measuring System

Test weight - nominal value	100 g
Test weight - volume (cm ³)	12,5
Test weight - volume uncertainty (cm ³)	0,1563
Permissible tolerance according to OIML (+/-)	0,5 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

Construction

Shape weight	Knob
Dimensions (Ø×H)	22×39 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	22×22×39 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg