

KERN 327-08

KERN

Individual weight, OIML Class F1, 200 g, Knob, stainless steel polished (OIML)



Services (optional)

Article number for conformity assessment	952-438
Article number for DAkkS calibration	962-438
Article number for guaranteed plus-deviation	921-438

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

Measuring System

Test weight - nominal value	200 g
Test weight - volume (cm ³)	25
Test weight - volume uncertainty (cm ³)	0,3125
Permissible tolerance according to OIML (+/-)	1 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

Construction

Shape weight	Knob
Dimensions (Ø×H)	28×48,5 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	28×28×48 mm
Net weight	0,2 kg
Shipping method	Parcel service
Net weight approx.	0,20 kg
Gross weight approx.	0,20 kg
Shipping weight	0,2 kg