

KERN 327-09

KERN

Individual weight, OIML Class F1, 500 g, Knob, stainless steel polished (OIML)



Services (optional)

Article number for conformity assessment	952-439
Article number for DAkkS calibration	962-439
Article number for guaranteed plus-deviation	921-439

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

Measuring System

Test weight - nominal value	500 g
Test weight - volume (cm ³)	62,5
Test weight - volume uncertainty (cm ³)	0,7813
Permissible tolerance according to OIML (+/-)	2,5 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

Construction

Shape weight	Knob
Dimensions (Ø×H)	38×66 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	38×38×66 mm
Net weight	0,5 kg
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,50 kg
Shipping weight	0,5 kg